


Part Number: **C30PU-RFV1872VAB**



*not to scale*

**Full Description:** Catheter for rat fem vein, PU 3Fr 25cm, collars @ 7.2, 12.2cm, VABR1B/22 attached

**Unit:** each      **Pouch Desc:** VAB FVC 3Fr L25 a7.2cm      **Pouch:** 4x10

**Approx Dead Vol:** 85  $\mu$ L      **Connects to:** 22ga      **EO Sterilized:**

**Intended Species:** rat      **Vessel:** femoral vein

**Usage Notes:** bead locations +/-1mm, length +/-3mm; volume estimate +/-15%

## TUBING

**Tip:** rounded      *Length 1*      *Length 2*  
 L1: 25 cm      L2:      cm  
**Material (IDxOD):** PU 3Fr (.62x1.02±.03mm)      **Material (IDxOD):**

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 7.2 cm	<input checked="" type="checkbox"/>
<b>Feature B:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 12.2 cm	<input type="checkbox"/>
<b>Feature C:</b>		c:      cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>	Rat VAB, 1ch, 22ga (VABR1B/22)	z:      cm	<input checked="" type="checkbox"/>
<b>Feature Y:</b>		y:      cm	<input type="checkbox"/>
<b>Feature X:</b>		x:      cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:** PNP3M injector      **Part 2:**