


Part Number: **C30PU-RIP1619VAB**



not to scale

Full Description: Catheter for rat peritoneum, PU 3Fr 20cm, perf holes, collars @ 0.7, 1.5, 4.0cm, VAB95BS attached

Unit: each **Pouch Desc:** VAB w/IPC PU 3Fr 20cm **Pouch:** 4x10

Approx Dead Vol: 72 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** peritoneum

Usage Notes: catheter UV bonded to VAB

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 20 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end

Feature	Description	Location	Fixed
Feature A:	Perfusion holes, 2 sets, 2mm apart	a: 0.2 cm	<input checked="" type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 0.7 cm	<input type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 1.5 cm	<input type="checkbox"/>

Proximal end

Feature Z:	Rat VAB, 22ga connector	z: 0.0 cm	<input checked="" type="checkbox"/>
Feature Y:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	y: 16.0 cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: VAH6M injector **Part 2:**