


Part Number: **C30PU-RIP1621VAB**



not to scale

Full Description: Catheter for rat peritoneum, PU 3Fr 30cm, perf holes, collars @ 0.7, 1.5, 4.0cm, VAB95BS attached

Unit: each **Pouch Desc:** VAB w/IPC PU 3Fr 30cm **Pouch:** 4x10

Approx Dead Vol: 105 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** peritoneum

Usage Notes: Catheter UV bonded to VAB. NOTE: CRL does not recommend 30cm catheter - too long.

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 30 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

| Distal (tip) end | | Location | Fixed |
|------------------|--|------------|-------------------------------------|
| Feature A: | Perfusion holes, 2 sets, 2mm apart | a: 0.2 cm | <input checked="" type="checkbox"/> |
| Feature B: | Small silicone collar for 3Fr (1.3mm OD), 2mm wide | b: 0.7 cm | <input type="checkbox"/> |
| Feature C: | Small silicone collar for 3Fr (1.3mm OD), 2mm wide | c: 1.5 cm | <input type="checkbox"/> |
| Proximal end | | | |
| Feature Z: | Rat VAB, 22ga connector | z: 0.0 cm | <input checked="" type="checkbox"/> |
| Feature Y: | Small silicone collar for 3Fr (1.3mm OD), 2mm wide | y: 16.0 cm | <input type="checkbox"/> |
| Feature X: | | x: cm | <input type="checkbox"/> |

PARTS ADDED TO POUCH

Part 1: VAH6M injector **Part 2:** cm