


Part Number: **C30PU-RJV1305**



*not to scale*

**Full Description:** Catheter for rat jugular vein, PU 3Fr 12cm, collar @ 3.0cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat JVC PU 3Fr 12cm **Pouch:** 3x7.5

**Approx Dead Vol:** 35  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** jugular vein

**Usage Notes:**

## TUBING

**Tip:** rounded *Length 1* *Length 2*  
**L1:** 12 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):**

## FEATURES

*Distal (tip) end* *Location* *Fixed*

**Feature A:** Small silicone collar for 3Fr (1.3mm OD), 2mm wide **a:** 2.5 cm

**Feature B:** **b:** cm

**Feature C:** **c:** cm

*Proximal end*

**Feature Z:** **z:** cm

**Feature Y:** **y:** cm

**Feature X:** **x:** cm

## PARTS ADDED TO POUCH

**Part 1:** **Part 2:**