


Part Number: **C30PU-RJV1445**



not to scale

Full Description: Catheter for rat jugular vein, PU 3Fr 12cm, round tip, fixed and moveable collar, perfusion holes, sleeve. Fits 22ga.

Unit: each **Pouch Desc:** Rat JVC 3Fr PU 12cm **Pouch:** 4x10

Approx Dead Vol: 40 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** jugular vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 12 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Perfusion holes, 2 sets, 2mm apart	a: 0.2 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 1.0 cm	<input checked="" type="checkbox"/>
Feature C:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	c: 2.0 cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:	Silicone sleeve for 3Fr (2.2mm OD), 15mm long	z: 1.0 cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**