


Part Number: **C30PU-RJV1462P**



*not to scale*

**Full Description:** Catheter for rat jugular vein, PU 3Fr 15cm, collars @ 2.5, 3.0cm, 2 PE dumbbells, 22ga white PinPort

**Unit:** each **Pouch Desc:** Rat JVC 3Fr PU 15cm PP white **Pouch:** 4x10

**Approx Dead Vol:** 50  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** jugular vein

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 15 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):**

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	<b>a:</b> 2.5 cm	<input type="checkbox"/>
<b>Feature B:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	<b>b:</b> 3.0 cm	<input type="checkbox"/>
<b>Feature C:</b>		<b>c:</b> cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>	PinPort, white, 22ga connector	<b>z:</b> 0.0 cm	<input type="checkbox"/>
<b>Feature Y:</b>	PE spool collar for 3Fr (1.6mm OD), 4 mm wide	<b>y:</b> 2.0 cm	<input type="checkbox"/>
<b>Feature X:</b>	PE spool collar for 3Fr (1.6mm OD), 4 mm wide	<b>x:</b> 3 cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:** PNP3M PinPort injector **Part 2:**