


Part Number: **C30PU-RJV1465**



not to scale

Full Description: Catheter for rat jugular vein, PU 3Fr 20cm, collars @ 2.5, 2.8cm, 19cm. Fits 22ga.

Unit: each **Pouch Desc:** Rat JVC PU 3Fr 20cm **Pouch:** 3x7.5

Approx Dead Vol: 60 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** jugular vein

Usage Notes: same as C30PU-RJV1405 but with collar at proximal end for identification

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 20 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.5 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 2.8 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	z: 1.0 cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**