


Part Number: **C30PU-RJV1508P**



not to scale

Full Description: Catheter for rat jugular vein, PU 3Fr 15cm, collars @ 2.5, 3.0cm, 2 PE dumbbells, 22ga red PinPort

Unit: each **Pouch Desc:** Rat JVC 3Fr PU 15cm PP red **Pouch:** 4x10

Approx Dead Vol: 50 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** jugular vein

Usage Notes:

TUBING

Tip: rounded **Length 1** **Length 2**
L1: 15 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.5 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 3.0 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
Feature Z:	PinPort, red, 22ga connector	z: 0.0 cm	<input type="checkbox"/>
Feature Y:	PE spool collar for 3Fr (1.6mm OD), 4 mm wide	y: 2.0 cm	<input type="checkbox"/>
Feature X:	PE spool collar for 3Fr (1.6mm OD), 4 mm wide	x: 3 cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: PNP3M PinPort injector **Part 2:**