


Part Number: **C30PU-RJV1821**



*not to scale*

**Full Description:** Catheter for rat jugular vein, PU 3Fr 20cm, collars @ 2.5, 2.8cm, 19cm, sleeve. Fits 22ga.

**Unit:** each **Pouch Desc:** Rat JVC PU 3Fr 20cm **Pouch:** 3x7.5

**Approx Dead Vol:** 60  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** rat **Vessel:** jugular vein

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 20 cm **L2:** cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm) **Material (IDxOD):** cm

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.5 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 2.8 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Silicone sleeve for 3Fr (2.2mm OD), 15mm long	z: cm	<input type="checkbox"/>
Feature Y:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	y: 2.0 cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**