


Part Number: C30PU-RJV2323



not to scale

Full Description:

Catheter for rat jug vein, PU 3Fr 20cm, 2.3cm insertion length (a2.1), rounded tip. 22ga luer stub.

Unit:

each

Pouch Desc:

Rat JVC 3Fr PU 20cm a2.1

Pouch:

4x10

Approx Dead Vol:

60

μL

Connects to:

22ga

EO Sterilized:

☒

Intended Species:

rat

Vessel:

jugular vein

Usage Notes:

For 225g SD rats when following Instech/A Lapierre surgical SOP (insertion point 5mm from clavicle)

TUBING

Tip:

rounded

Length 1

L1:

20

cm

Length 2

L2:

cm

Material (IDxOD):

PU 3Fr (.62x1.02±.03mm)

Material (IDxOD):

FEATURES

Distal (tip) end

Feature A:

Small silicone collar for 3Fr (1.3mm OD), 2mm wide

Location

a:

2.1

cm

Fixed

☐

Feature B:

b:

cm

☐

Feature C:

c:

cm

☐

Proximal end

Feature Z:

z:

cm

☐

Feature Y:

y:

cm

☐

Feature X:

x:

cm

☐

PARTS ADDED TO POUCH

Part 1:

LS22

Part 2: