


Part Number: C30PU-RJV2329



not to scale

Full Description: Catheter for rat jug vein, PU 3Fr 20cm, 2.9cm insertion length (a2.7), rounded tip. 22ga luer stub.

Unit: each

Pouch Desc: Rat JVC 3Fr PU 20cm a2.7

Pouch: 4x10

Approx Dead Vol: 60 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: rat

Vessel: jugular vein

Usage Notes: For 300g SD rats when following Instech/A Lapierre surgical SOP (insertion point 5mm from clavicle)

TUBING

Tip: rounded

Length 1

L1: 20 cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)

Length 2

L2: cm

Material (IDxOD):

FEATURES

Distal (tip) end

	Location	Fixed
Feature A: Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 2.7 cm	☐
Feature B:	b: cm	☐
Feature C:	c: cm	☐

Proximal end

Feature Z:	z: cm	☐
Feature Y:	y: cm	☐
Feature X:	x: cm	☐

PARTS ADDED TO POUCH

Part 1: LS22

Part 2: