


Part Number: C30PU-RJV2332



not to scale

Full Description:
Catheter for rat jug vein, PU 3Fr 20cm, 3.2cm insertion length (a3.0), rounded tip. 22ga luer stub.

Unit: each

Pouch Desc: Rat JVC 3Fr PU 20cm a3.0

Pouch: 4x10

Approx Dead Vol: 60 μ L

Connects to: 22ga

EO Sterilized: ☒

Intended Species: rat

Vessel: jugular vein

Usage Notes:
For 350g SD rats when following Instech/A Lapierre surgical SOP (insertion point 5mm from clavicle)

TUBING

Tip: rounded

Length 1
L1: 20 cm

Length 2
L2: cm

Material (IDxOD): PU 3Fr (.62x1.02 $\pm$.03mm)
Material (IDxOD):

FEATURES

Distal (tip) end

Feature A: Small silicone collar for 3Fr (1.3mm OD), 2mm wide

Feature B:

Feature C:

Location
a: 3.0 cm
b: cm
c: cm

Fixed
☐
☐
☐

Proximal end

Feature Z:

Feature Y:

Feature X:

z: cm
y: cm
x: cm

☐
☐
☐

PARTS ADDED TO POUCH

Part 1: LS22
Part 2: