

CM375KRP

Snap-in spring swivel arm for rat shoebox cage top and plastic swivel



SPECIFICATIONS

Compatible swivel diameter	.375in, 9.5mm (Instech 375/22PS)
Compatible swivel weight	2-3g (plastic, not stainless steel)
Rod height	2.5in, 6.4cm
Spring arm length	5.5in, 14cm
Clearance required above cage lid	7.5in, 19cm
Compatible cage lid wire spacing	9-11mm (7-9mm openings)
Cage wash compatible	Yes (remove swivel from clip)
Assembly tools required	Phillips head screwdriver
Intended use	Tethered rat infusion or sampling (spring forces not appropriate for mice)

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/mounts/snap-in-spring>

Instructions For Use

Product Video <https://www.instechlabs.com/how-to-mount-cm375krp>

NOTE: All products are designed and sold for research only; not for human use.
