

# CMMI/65

Mouse micro-isolator cage mount



## SPECIFICATIONS

Boom arm length	10cm (pivot to center of swivel)
Height above cage edge	6cm for bracket, 10cm with arm upright
Compatible swivels	Instech plastic or stainless steel, 0.375in diameter
Swivel mount type	Gimbal with slot and 0.050in set screw
Counter-balance range	0-40g with adjustable spring tension
Range of vertical motion	Approximately 12cm
Horizontal range of motion	Can rotate freely or be fixed with .050in set screw
Vertical collar adjustment range	Approximately 3cm
Intended use	Holding Instech fluid swivels inside mouse microisolator cages

## FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/mounts/micro-isolator-cage>

Instructions For Use

Product Video

NOTE: All products are designed and sold for research only; not for human use.

---