

MCS/5A

Five Channel Swivel: 4 Microdialysis, 1 22ga Channels



SPECIFICATIONS

Number of channels	5
Dead volume - microdialysis channels	6-7µL
Dead volume - 22ga channels	15-20µL
ID - microdialysis channels	0.009in (0.2mm)
ID - 22ga channels	0.015in (0.4mm)
OD - all channels	0.028in (0.7mm) - 22ga
Materials in fluid path	Titanium, polyamide, Teflon
Rotational activity output	Analog signal, ±2V max
Lever arm length	5.75in (14.6cm)
Swivel dimensions (LxD)	7x0.6in (18x1.5cm)
Swivel weight	95g, counterbalanced

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/swivels-mounts/swivels/five-channel>

Instructions For Use https://www.instechlabs.com/hubfs/pdfs/manuals/7700050178_MCSmanual.pdf

Product Video

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA
www.instechlabs.com