

P400/MS

Mounting system for Instech model 400 syringe pump



SPECIFICATIONS

Contents	C-clamp with V-lock base station, aluminum pump mounting block with rubber feet, V-lock wedge, (2) 1/4-20 pump mounting screws, (3) 4-40 3/8in screws for V-lock wedge, screw driver, 5/32 allen key
Attaches to	Shelf, bench, pole
C-clamp maximum opening	42mm
Orientation	Vertical or horizontal, multiple configurations possible
Materials	Aluminum, plastic, rubber
Weight, installed	0.35kg
Package size	9x4x3in
Packaged weight	0.5kg (1lb)
Intended use	Mounting of Instech P400 syringe pumps

FOR MORE INFORMATION

Webpage <https://www.instechlabs.com/products/pumps/instech-syringe-pumps/instech-syringe-pump/pump-mounts>

Instructions For Use https://www.instechlabs.com/hubfs/pdfs/manuals/p400-ms_ifu.pdf

Product Video

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA
www.instechlabs.com