

PNP3M

PinPort injectors, 50 sterile pouches of 5



SPECIFICATIONS

Dead volume	~30 μ L, mostly in luer gap
Size of tube puncturing PinPort	25ga (.010in, .25mm ID)
Materials in fluid path	Stainless steel, polycarbonate
Sterility	Sold sterilized by ETO gas; single use only; do not autoclave
Compatible syringes	Any with male slip or locking luer connector compliant with ISO 80369-7
Compatible ports	For use with Instech PinPorts only

FOR MORE INFORMATION

Webpage </products/tubing-connectors-pinports/pinports/replacement-injectors>

Instructions For Use

Product Video

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com