

VAHR1H/22-1P

Rat VAH harness, 22ga PinPort, with side-mounted red 22ga PinPort



SPECIFICATIONS

Port configuration	White port for tether, side mounted red port for intermittent access
Port diameter	3mm (PinPort™)
Port dead volume	3µL
Catheter connectors	22ga (.028in, .71mm OD)
Compatible catheters	3Fr (.60-.64mm ID) polyurethane
Mating tethers	VAHR1T/22, KVAHR1T/22 or similar (not compatible with VAH95T designs)
Mating injector	PNP3M (not compatible with VAH6M)
Materials in fluid path	Acetal, stainless steel, silicone
Saddle dimensions	1.13x1.13in (2.9cm)
Belly band length	9in (23cm)
Weight	6g
Typical application	Rat infusion (continuous) and blood sampling (intermittent)

FOR MORE INFORMATION

[Webpage](#)

[Instructions For Use](#)

[Product Video](#)

NOTE: All products are designed and sold for research only; not for human use.

Instech Laboratories, Inc.

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

www.instechlabs.com