

## VAHR1H/22

Rat VAH harness, 22ga PinPort



### SPECIFICATIONS

<b>Part diameter</b>	3mm (PinPort™)
<b>Port dead volume</b>	3µL
<b>Catheter connector</b>	22ga (.028in, .71mm OD)
<b>Compatible catheters</b>	3Fr (.60-.64mm ID) polyurethane
<b>Mating tethers</b>	VAHR1T/22, KVAHR1T/22 or similar (not compatible with VAH95T designs)
<b>Mating injector</b>	PNP3M (not compatible with VAH6M)
<b>Materials in fluid path</b>	Acetal, stainless steel, silicone
<b>Saddle dimensions</b>	1.13x1.13in (2.9cm)
<b>Belly band length</b>	9in (23cm)
<b>Weight</b>	5g
<b>Included in pouch with harness</b>	Blue port cover, PNP3M injector
<b>Typical application</b>	Rat infusion

### FOR MORE INFORMATION

[Webpage](#)

[Instructions For Use](#)

[Product Video](#)

NOTE: All products are designed and sold for research only; not for human use.

---

**Instech Laboratories, Inc.**

450 South Gravers Road | Plymouth Meeting, PA 19462 | USA

[www.instechlabs.com](http://www.instechlabs.com)