

Blood Sampling and Catheter Maintenance Kits

Part Number	Syringe No.	Heparin (USP units/mL)	Saline (NaCl Conc.)	Glycerol Conc.	Syringe Type	Fill Volume (mL)	PinPort™ Injector
USPK-S-1-HGS500-0.3-NJ1P	1st		0.9%		B-D® 1mL	1	✓
	2nd	500		50%	HSW NORM-JECT® 1mL	0.3	✓
	3rd				B-D® 1mL	0	✓
USPK-S-0.5-HS020-1-NJ1P	1st		0.9%		HSW NORM-JECT® 1mL	0.5	✓
	2nd	20	0.9%		HSW NORM-JECT® 1mL	1	✓
	3rd				HSW NORM-JECT® 1mL	0	✓


PRODUCT PAGE



Click or Scan

INSTRUCTIONS

FOR USE



Click or Scan

INSTECH

Instech Laboratories, Inc.

450 South Gravers Road, Plymouth Meeting, PA 19462 | USA

www.instechlabs.com