

# INSTECH STAINLESS STEEL SWIVEL REPAIR KITS

REPAIR KIT PART NUMBER	375R/D/xx	375R/xx	375R/D/TK	375R/CSxx	375R/D/SSxx	375R/BR
OLD PART NUMBER	125/D/xx		125/D/TOOLKIT	125/xx	SIDSEAL/xx	125/BEARINGS
DESCRIPTION	Tools and parts to repair one dual channel swivel (xx=gauge)	Tools and parts to refurbish five single channel swivels (xx=gauge)	Tools for refurbishing dual channel swivels (no parts included)	Parts to repair center channel on 5 dual or single channel swivels	Parts to repair sideseals on five dual channel swivels	Parts to repair bearings on five dual channel swivels
REPAIR TOOLS INCLUDED						
Tweezers	1		1			
8ga tube (for removing bearings from outlet tube)	1	1	1			
12ga tube (for removing outlet tube from main body)	1		1			
Bearing removal base	1		1			
Polishing paper	1	1		1	1	1
Repair manual	1	1	1	1	1	1
REPLACEMENT PARTS INCLUDED						
Center channel seal (xx=size)	1	5		5		
Sideseal body with sideseals (xx=size)	1				5	
Ball bearings for dual (2 per swivel)	2					10
Torlon bearings for single (1 per swivel)		5				
1/4" clip rings (for main body)	3				10	
1/8" clip rings (for bearings)	3	15			15	15

For prices see: <http://www.instechlabs.com/Infusion/swivels/swivelrepair.php>

## KEY

For dual channel swivels (excluding QM, LT)
For single channel swivels
For singles or duals

# INSTECH

Instech Laboratories, Inc.  
Plymouth Meeting, PA USA  
TEL (800) 443-4227  
[www.instechlabs.com](http://www.instechlabs.com)