


Catheter Flush and Lock Prefilled Syringes

Part Number	Heparin (USP units/mL)	Dextrose Conc.	Glycerol Conc.	Saline (NaCl Conc.)	Syringe Type	Fill Volume (mL)	Syringe Qty.	Flush Solution	Lock Solution	PinPort™ Injector
USP-S-0.5-BD1P-50	0			0.9%	B-D® 1mL	0.5	50	✓	✓	✓
USP-S-1-BD1P-50	0			0.9%	B-D® 1mL	1	50	✓	✓	✓
USP-HS-005-0.5-NJ1P-50	5			0.9%	HSW NORM-JECT® 1mL	.5	50	✓	✓	✓
USP-HS-010-0.5-NJ1P-50	10			0.9%	HSW NORM-JECT® 1mL	0.5	50	✓	✓	✓
USP-HS-010-1-NJ1P-50	10			0.9%	HSW NORM-JECT® 1mL	1	50	✓	✓	✓
USP-HS-020-0.5-NJ1-50	20			0.9%	HSW NORM-JECT® 1mL	0.5	50	✓	✓	
USP-HS-020-0.5-NJ1P-50	20			0.9%	HSW NORM-JECT® 1mL	0.5	50	✓	✓	✓
USP-HS-020-1-NJ1P-50	20			0.9%	HSW NORM-JECT® 1mL	1	50	✓	✓	✓
USP-HS-250-1-NJ1P-50	250			0.9%	HSW NORM-JECT® 1mL	1	50		✓	✓
USP-HD25-100-0.5-NJ1P-50	100	25%			HSW NORM-JECT® 1mL	0.5	50		✓	✓


PRODUCT PAGE



Click or Scan

INSTRUCTIONS

FOR USE



Click or Scan

INSTECH

Instech Laboratories, Inc.

450 South Gravers Road, Plymouth Meeting, PA 19462 | USA

www.instechlabs.com